Colostrum 101

Everything you need to know about collecting, storing, and feeding your newborn their first nutrient-rich milk.

Why collect colostrum?

One of the most helpful reasons to collect colostrum is to avoid formula supplementation.* It can be especially beneficial for moms who:

- Have gestational diabetes
- Are carrying twins or multiples
- Have pregnancy-induced hypertension
- Have a baby who is not effectively latching
- Have a baby who spends time in a special care unit

*Always consult your provider before collecting colostrum

By collecting colostrum prenatally and freezing it, parents can ensure that, regardless of what happens in the first few days, there is plenty of “liquid gold” to supplement until breastfeeding can resume.

Colostrum has many benefits

- Helps protect your baby from illness
- Helps your baby’s developing digestive system learn how to process food
- Acts as a diuretic to remove meconium

What’s inside

1 Collection Cup

Comfortably and easily catch hand-expressed colostrum

20 Syringes (10x 1 ml and 10x 5 ml)

Pre-sterile & individually packaged to store small amounts

22 Labels

Easily organize and track each drop with labels

2 Storage Boxes

Safely store each syringe of colostrum in the refrigerator or freezer

PLUS! This colostrum education guide, created in partnership with pediatrician and breastfeeding medicine specialist, Dr. Laurie Jones, on the What, Why, When, and How of collecting, storing, and feeding.

Intended Use

The Colostrum Collection Kit is intended for collecting, storing, and freezing colostrum from 37 weeks of pregnancy. It enables feeding your newborn breast milk if your baby’s provider recommends supplementation.

Everyone should consult with their provider before using this product, especially if you:

- Have a cervical suture (stitches) in place
- Have experienced premature labor in the past or a history of miscarriage
- Have had contractions, vaginal bleeding and/or premature rupture of membranes (breaking early) during this pregnancy
- Have placenta previa or accreta

⚠️ WARNINGS

Using an electric or manual breast pump too early can increase the risk of premature labor.

Before each use, visually inspect parts. If damage is found, discontinue use until parts have been replaced.

To prevent contamination, always wash your hands prior to handling parts and disinfect the collection cup prior to every use. If you experience contractions that are painful, regular, or accompanied by cramping or bleeding, stop immediately and seek medical attention.

CHOKING HAZARD! Keep out of reach of children under the age of three.

How it works

Learn

Know when to collect colostrum, how to store it, and when to feed your baby

Collect

In the final weeks of pregnancy, start hand expressing into the collection cup

Store

Safely and effectively freeze in the storage box in 1ml or 5ml labeled syringes

Feed

Talk with your provider about how and when to feed your newborn

Read on for step-by-step directions.

What is colostrum?

Colostrum is the first, most concentrated form of breast milk, commonly referred to as “liquid gold.” It is rich in antibodies and immune factors that help protect your baby. Your body may begin producing colostrum as early as 12–18 weeks in pregnancy and will continue to make it until your milk transitions 3–5 days after birth.

When can I start collecting colostrum?

Anyone with an otherwise uncomplicated pregnancy can benefit from collecting colostrum at or after 37 weeks gestation and post-birth before your transitional milk comes in. Always talk to your provider before you begin and only hand-express colostrum.

How much and how often should I collect colostrum?

If you collect prenatally at 37 weeks, it’s recommended to practice once a day until your baby arrives. Do not express more than three times in 24 hours. Before delivery, you can expect to produce about 1.5–2 milliliters per session.

First-time moms may not get much, especially in the beginning. By design, your body will only produce a low volume to accommodate your newborn’s small stomach for the first several days.

How much colostrum will my baby need?

Colostrum packs a major nutritional punch for your baby, as a small amount will contain all the nourishment, fat, and immune factors your baby needs in the first several days of their life.

For example, in the first 24 hours of your baby’s life, they may receive about 15ml (about half an ounce) of colostrum.

How do I express colostrum?


Learning to express colostrum by hand is recommended by the American College of Obstetricians and Gynecologists. Research shows that it can help you save every last drop versus trying to pump into a container.

Step-by-step guide to expressing colostrum

Before collecting colostrum, boil the collection cup for five minutes and dry completely. If the syringe packaging is damaged, do not use.


Step 1
With clean hands, gently massage your breast.

Step 2
Next, position your thumb and index finger in the area between your nipple and areola. They should be in a straight line.

Step 3
Press back into the chest wall and simultaneously compress.


Step 4
Repeat a few times. If you are unable to express milk, move your fingertips closer to your nipple or further away to find the "sweet spot" where the milk is expressed. Catch colostrum in the collection cup.

Step 5
Relax and repeat until milk flow slows. Then rotate and adjust hands to a different area of the breast and begin again.


How to store colostrum

Step 1
Remove the syringe from its pre-sterile packaging. Take the cap off, but don't dispose of it; you'll need it for storage. Use the syringe to draw up colostrum from the collection cup. Do not fill past the max capacity line (1ml or 5ml). Remove all excess air from the syringe to ensure maximum retention of colostrum nutrients.

Step 2
Add the cap and label with your full name, birthdate, and date collected.

Step 3
Place the syringe in the storage box and freeze.

Ways to thaw colostrum

- For immediate use, thaw under running water and use within two hours. For next day use, thaw in the fridge for 12 hours and use within 24 hours
- Place the syringe in a bowl of warm water. Swirling the syringe in the water is an effective way to thaw it
- Run warm (not hot) water over it
- If you use a bottle warmer, ensure it has a defrost setting and is not on a high or steam setting, as this can damage the milk or overheat it

Storage guidelines

Colostrum can be stored at room temperature for up to four hours or in the refrigerator for up to four days. You can store colostrum in the freezer for 12 months but this liquid gold is most useful for baby within the first six months of life.

IMPORTANT!

- Never microwave to thaw or heat colostrum (breast milk). It can create hot spots that can burn your baby and destroy nutrients
- Use thawed milk within 24 hours. If it is left out at room temperature, it must be used within two hours
- Once thawed, milk cannot be frozen again
- If your baby doesn't finish a feeding, discard the remaining milk within 1–2 hours
- When in doubt, throw it out: If the milk appears stringy or foul, discard it

When and how to feed colostrum

At the hospital

If supplementation is needed for any reason during your hospital stay and you prefer not to use formula, talk to your provider about their preferred method of feeding colostrum to infants. Common methods include feeding by spoon, syringe, bottle, and cup.

If your provider deems supplementation is necessary

If you have a delay of copious milk production, you and your provider may decide to supplement with colostrum.

As an antibody booster

If you choose not to feed your baby colostrum in the first few days postpartum, you can add it to purées or a bottle to boost nutrition or apply it to irritated skin as a natural healing method. It is recommended to use it in the first six months.

Remember

If you are supplementing, be sure that any colostrum you collect goes to the baby. Don't worry about boosting your supply to the point of saving extra.

In the first few days and weeks, just focus on ensuring your baby is getting all the milk you are making.

The amount of colostrum collected prenatally does not indicate future milk supply.

Finally, don't feel pressured. Collection is optional and not necessary for all parents.

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